EV9627-14 2SR-BW Filter Cartridge		
Controls Scale Formation in RO Applications		
<text><text><text><text></text></text></text></text>	2SR-BW Cartridge: EV9627-14	BENEFITS For exclusive use as pre- treatment to an RO Membrane. Controls scale formation to extend membrane life. Sanitary cartridge replacement is simple, quick and clean. Internal filter parts are never exposed to handling or contamination.
INSTALLATION TIPS Never use saddle valve for connection Use 3/8" water line Install vertically with cartridge hanging down Allow 2-1/2" clearance below the cartridge for easy cartridge replacement Flush 2SR-BW cartridge by running water through filter for three minutes at full flow	OPERATION TIPS Change cartridges on a regular 6 month preventative maintenance program Change cartridges when capacity is reached or when pressure falls below 10 psi Service flow rate must not exceed 1.0 gpm Always flush the filter cartridge at time of installation and cartridge change	A P P L I C A T I O N / S I Z I N G For reverse osmosis applications Rated Capacity: 25,000 gallons (94,635l)

2SR-BW Filter Cartridge

S P E C I F I C A T I O N S

Overall Dimensions: 10.25"H x 3.25"D

Service Flow Rate: Maximum 1.0 gpm (3.8 lpm)


Rated Capacity: 25,000 gallons (94,635l)

Pressure Requirements: 10 - 125 psi (0.7 - 8.6 bar), non-shock

Maximum Temperature: 35 - 100°F (2 - 38°C)

No electrical connection required

The contaminants or other substances removed or reduced by this drinking water system are not necessarily in your water. Do not use with water that is microbiologically unsafe or of unknown quality without adequate disinfection before or after the system.



WARRANTY

Everpure replaceable elements (filter cartridges and water treatment cartridges) are covered by a limited warranty against defects in material and workmanship for a period of one year after date of purchase. See printed warranty for details. Everpure will provide a copy of the warranty upon request.





Liberty Vertriebs GmbH Bonner Straße 103 40589 Düsseldorf

H Tel.: 0211-753520 Fax.: 0211-753575

Internet: www.libertyshop.de Email: verkauf@libertyshop.de